

DIETARY SUPPLEMENT

Helps Calm and Soothe Babies\*

# Sucro 24

STERILE SUCROSE AT 24% FROM ACER SACCHARUM SAP

\*This statement has not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.

100 X 2 ML (0.066 fl.oz.)

### Supplement Facts

20 single serve units per container - 100 containers per pack  
Serving size 2 drops (0.003 fl.oz.)

Amount per serving  
Total Carbohydrate Less than 0.05 g from Acer saccharum Sap

Not a significant source of total sugars

OTHER INGREDIENT: Purified water.

#### SUGGESTED USE (INFANT 0-12 MONTHS):

■ Use according to hospital protocol. ■ 2-3 servings 1 to 10 times a day. ■ Sucro24 works better when administered at one (1) to two (2) minutes prior to the discomforting procedure. ■ Sucro24 can be used even if the baby isn't able to eat or drink anything before a procedure. ■ Sucro24 will not affect blood test results.

#### DIRECTIONS FOR USE:

Pour two (2) drops (0.003 fl.oz, in total) at the tip of the baby's tongue. Administer two (2) servings two (2) minutes apart. A third (3rd) serving can be administered after a two (2) minute delay between each serving.

#### CAUTIONS AND WARNINGS:

■ NOT INJECTABLE. ■ Do not use if seal is broken or if the bottle is damaged. ■ Ready-to-use solution. Do not dilute. ■ Consult a healthcare practitioner prior to use in case of gastrointestinal disorder (such as necrotizing enterocolitis), hemodynamic instability, esophageal artesia or tracheoesophageal fistula, severe neurological distress (severe asphyxia), severe respiratory distress of unstable sugar level. ■ Do not use on premature or low birth weight newborns without permission from a pediatrician.

DIETARY SUPPLEMENT

Helps Calm and Soothe Babies\*

# Sucro 24

STERILE SUCROSE AT 24% FROM ACER SACCHARUM SAP

\*This statement has not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.

100 X 2 ML (0.066 fl.oz.)

# Sucro 24

100 X 1 ML (0.033 fl.oz)

inopro 

MADE IN CANADA  
by Inopro inc.  
Cowansville, QC  
Canada J2K 2Y3  
www.inopro.ca  
514 434-8481



Storage: 40° - 80°F  
Code: PFC-1002USA.1

Lot #   
Best By Date: 

