

**SUGGESTED USE:** Add 7 drops per pint of milk, refrigerate for 24 hours, then use as desired. If lactose intolerance symptoms persist, refrigerate 48 hours before using or add more drops. Alternatively, add 7 drops directly to milk or any dairy product just prior to drinking or eating, or use as directed by your healthcare professional.

Do not use if tamper seal is damaged.  
Keep out of reach of children.

**STORAGE:** Refrigerate after opening.  
Best consumed within 60 days after opening.

**DOES NOT CONTAIN:** Milk, eggs, fish, shellfish, tree nuts, peanuts, wheat, soy, gluten, or GMO, artificial colors, artificial flavors, or artificial preservatives.

**Lactase Drops Meets or Exceeds cGMP Quality Standards.**

Lactase Drops - 1/2 FL OZ - NEW



X000U9501R

TARGETED SUPPORT®



SeekingHealth®

# Lactase Drops

Liquid Lactase  
Enzyme Supplement

1/2 FL OZ (15 ML)

## Supplement Facts

Serving Size: 7 Drops

Servings Per Container: 52

	Amount Per Serving	%DV
Lactase	750 Lactase Units	*

\*Daily Value (DV) not established.

Other Ingredients: Glycerin and purified water.

**Manufactured for and Distributed by:**

Seeking Health

3140 Mercer Ave

Bellingham, WA 98225

(800) 547-9812

Find Other Effective Health Products At

[www.seekinghealth.com](http://www.seekinghealth.com)

LACT-HALFOZ-LIQ-A REV. 022422

**MANUFACTURED IN THE USA**  
FROM GLOBALLY SOURCED INGREDIENTS  
THIRD-PARTY INDEPENDENTLY TESTED