

DIRECTIONS: Take one softgel daily, or use as directed by your healthcare professional.

Consult your healthcare professional before use. Individuals taking blood thinners or other medication should discuss potential interactions with their healthcare professional. Do not use if tamper seal is damaged.

OmegaPure PRM features highly concentrated fish oil with a standardized quantity of the pro-resolving mediators (PRMs) 18-HEPE, 17-HDHA, and 14-HDHA. Metabolites of PRMs, known as specialized PRMs, provide building blocks to support the natural resolution of the immune response.*

*These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.

STORAGE: Keep closed in a cool, dry place out of reach of children.

FORMULATED TO EXCLUDE: Wheat, gluten, corn, yeast, soy, dairy products, crustacean shellfish, peanuts, tree nuts, egg, sesame, ingredients derived from genetically modified organisms (GMOs), artificial colors, and artificial sweeteners.

ePV:



200092485
© 2024 XYMOGEN
REV. 110824



OmegaPure PRM

Featuring Pro-Resolving Mediators*

PROFESSIONAL FORMULAS

60 Softgels | Dietary Supplement

Supplement Facts

Serving Size: 1 Softgel
Servings Per Container: 60

Amount Per Serving %Daily Value		
Calories	10	
Total Fat	1 g	1%†
Cholesterol	5 mg	2%
Fish Oil Concentrate	800 mg	**
Total Omega-3 Fatty Acids	590 mg	**
EPA (eicosapentaenoic acid)	245 mg	**
DHA (docosahexaenoic acid)	195 mg	**
Total Pro-Resolving Mediators	400 mcg	**
18-HEPE (18-hydroxyeicosapentaenoic acid)		
17-HDHA (17-hydroxydocosahexaenoic acid)		
14-HDHA (14-hydroxydocosahexaenoic acid)		

† Percent Daily Values are based on a 2,000 calorie diet.

** Daily Value not established.

Other Ingredients: Softgel (bovine gelatin, glycerin, purified water), and natural lemon flavor.

Contains: Fish (anchovy, sardine, herring, and mackerel) and squid.

Distributed by:
XYMOGEN®
Professional Formulas
6900 Kingspointe Pkwy. Orlando, FL 32819
xymogen.com 800-647-6100



Scan to learn